
T.A.C. PIN®

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
• Color code for size identification:
 Thread length 15 mm
 Thread length 20 mm
 Thread length 25 mm

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• Lot number and reference marking for easier product traceability.

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Titanium pin with compression thread.

- Variable thread pitch for true and adapted compression.
- Direct bone implantation.
 - NB :** According to the bone quality, a pre-drill may be required.
- 3 length of threads : 15, 20 and 25 mm.
- Material in titanium alloy TiAl6V4 ISO 5832-3 ASTM F136.
- Patented implant.

IMPLANT SELECTION

The choice of the implant is depending from the osteosynthesis to perform. The ideal choice is the implant which offers most threads at both sides from the osteotomy or fracture with the short non-threaded lag part at least at the level of the osteotomy / fracture. The diameter at the non-threaded part is 1,2 mm.

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OTHER POSSIBLE INDICATIONS

- Arthrodesis in hand surgery.
- Mono or bi-cortical osteotomies in the hand.
- Distal or proximal metacarpal osteotomies.

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The tools needed for setting a T.A.C.'pin® are :

- a wire driver accepting a 1,6mm diameter pin
- an optional universal pin cutter.

The pin should be inserted with the wire driver in a speed as slow as possible, as corresponding to manual drilling. High speed or / and fixed speed drills are not advised.


No angular force should be applied on the pin in order to prevent mal positioning of the same of fragments or bending of the T.A.C.'pin®.

According to the bone quality, a drilling of the cortex with a 1mm K-wire may be required for easier setting of the T.A.C.'pin®.

It is recommended to control and assist the ideal location and the formal setting of the T.A.C.'pin® by fluoroscopy or X-rays. The non-threaded part should be ideally at the level of the fracture / osteotomy. It may pass this level, but the whole threads should not be only in the distal bone fragment.

The pin may be cut off with a sharp wire cutter without bending the remaining implanted part of the pin. A percutaneous T.A.C.'pin® would indeed be removed with a wire driver only as it is guaranteed to be straight. Removal of a bent pin may cause dislocation of the bone structures or fracture of the pin itself.

REFERENCE	T.A.C.' PIN® DIAM 1.6MM	DESCRIPTION
120 015	• LONG. FILETAGE 15MM • THREAD LENGTH 15MM	
120 020	• LONG. FILETAGE 20MM • THREAD LENGTH 20MM	
120 025	• LONG. FILETAGE 25MM • THREAD LENGTH 25MM	

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THE T.A.C.'PIN® IS TO BE IMPLANTED THANKS TO THE USE OF STANDARD SURGICAL POWER TOOL EQUIPMENT

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- The products are manufactured and referenced within the frame of the standards in force.
- Implantation procedures are described in the surgical technique.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- WARNING : Federal law (USA) restricts this device to sale by or on the order of a physician.

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NEWDEAL SAS