
STIMULAN®

Rapid Cure

Instructions For Use

Number I.025
Rev 062122

Ref Code	Description	Size	UDI-DI Number
620-003*	Stimulan Rapid Cure	3cc	15060155711451
620-005	Stimulan Rapid Cure	5cc	15060155711024
620-010	Stimulan Rapid Cure	10cc	15060155711031
620-020	Stimulan Rapid Cure	20cc	15060155711048

* CE Pending. Not yet available in EU, AUS and CAN

Sample Type	Average Days to Complete Dissolution ±SD (n=10)
Stimulan Rapid Cure	54 ± 8.3
Stimulan Rapid Cure Vancomycin	32.4 ± 2.8
Stimulan Rapid Cure Tobramycin	32 ± 3.0
Stimulan Rapid Cure Gentamicin	35.4 ± 2.6

Contents

(US) United States of America.....	2
(EN) English.....	4
(BA) Bosanski.....	6
(BG) български.....	8
(BR) Português.....	10
(CZ) Česky.....	12
(DK) Dansk.....	14
(DE) Deutsch.....	16
(EE) Eesti.....	18
(ES) Español.....	20
(FI) Suomi.....	22
(FR)Français.....	24
(GR) Ελληνικά.....	26
(HR) Hrvatski.....	28
(HU) Magyar.....	30
(IT) Italiano.....	32
(LT) Lietuva.....	34
(ME) Crnogorski.....	36
(MK) Македонски.....	38
(NL) Nederlands.....	40
(NO) Norsk.....	42
(PL) Polskie.....	44
(PT) Português.....	46
(RO)Română.....	48
(RS) Serbski.....	50
(SE) Svenska.....	52
(SI) Slovenski.....	54
(SK) Slovensky.....	56

(US) UNITED STATES ONLY

STIMULAN® for use in the presence of infection

PRODUCT DESCRIPTION AND MATERIALS

STIMULAN® Rapid Cure is provided sterile for single patient use. STIMULAN® Rapid Cure contains calcium sulfate powder and mixing solution in pre-measured quantities so that when mixed together in a sterile mixing bowl, the resultant paste is to be digitally packed into open bone void/gap to set in situ or placed into the mould provided, the mixture sets to form beads. The biodegradable, radiopaque beads are resorbed in approximately 30 – 60 days when used in accordance with the device labelling. STIMULAN® Rapid Cure is manufactured from synthetic implant grade calcium sulfate dihydrate ($\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$) that resorbs and is replaced with bone during the healing process. Also, as the bone void filler beads are biodegradable and biocompatible, they may be used at an infected site. Not made with natural rubber latex.

INTENDED USES

To fill bone void or defect.

INDICATIONS FOR USE

STIMULAN® Rapid Cure is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. These osseous defects may be surgically created or created from traumatic injury to the bone. STIMULAN® Rapid Cure is intended to be gently packed into bony voids or gaps of the skeletal system (i.e., extremities, pelvis, and posterolateral spine). STIMULAN® Rapid Cure provides a bone graft substitute that resorbs and is replaced with bone during the healing process. STIMULAN® Rapid Cure is biodegradable and biocompatible and may be used at an infected site.

To fill a bone void or defect created by:

- surgery
- a cyst
- a tumour
- osteomyelitis
- traumatic injury

CONTRAINDICATIONS

- filling of defects which are intrinsic to the stability of the bony structure
- severe vascular or neurological disease
- uncontrolled diabetes
- severe degenerative bone disease
- pregnancy
- uncooperative patient who can't or won't follow post-operative instructions including individuals who abuse drugs or alcohol.
- hypercalcaemia
- renal compromised patients

WARNINGS AND PRECAUTIONS

The device is for single use only and must not be reused. Do not resterilise the device. Do not use beyond expiration date on the label. Do not use if product package shows signs of tampering or damage. Do not add other substances to the device. Adding other substances may alter the safety and effectiveness of this product. The device must be implanted where surrounding bone is healthy and vascular but not actively bleeding at time of insertion. Do not use for bony voids in the cranium, due to the potential risk of contact with cerebrospinal fluid (CSF) or dura mater. Do not irrigate the site following implantation. The device is intended for manual application and is not intended for injection. Injection of the device should not be conducted as it could result in over-pressurisation, which may lead to device extrusion beyond the intended application site or to embolisation of fat or the device into the bloodstream. Avoid overfilling the bone void or pressurising the treatment site. Incomplete or inadequate soft tissue coverage may result in device migration and/or effusion or serous exudate. Do not implant unless adequate tissue coverage and/or containment can be achieved. The patient should be advised to report any related pain, swelling, fever or unusual incidences. The patient is to be cautioned to govern activities and protect the surgical site from unreasonable stresses and follow the instructions of the physician with respect to follow-up care and treatment. The patient is to be warned of surgical risks and made aware of possible adverse effects. Concurrent use of locally administered antibiotics may affect the setting time.

ADVERSE EFFECTS

Peripheral neuropathies have been reported following surgery. Sub-clinical nerve damage occurs more frequently, possibly the result of surgical trauma. Material sensitivity/allergic reactions in patients following surgery have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. Infection can lead to failure/removal of the void filler. Adverse effects may include but are not limited to: wound complications including hematoma, site drainage, bone fracture, infection and other complications that are possible with any surgery, fracture or extrusion of the bone void filler, with or without particulate debris generation, deformity of the bone at the site and incomplete or lack of osseous ingrowth into bone void as is possible with any bone void filler.

STERILISATION

When sterilisation dot is red product has been gamma irradiated.

Unless opened or damaged, the product is supplied sterile within a double pack. Check all packaging for punctures or other damage. Sterilisation of components is by gamma radiation to achieve a Sterility Assurance Level (SAL) of 10^{-6} . When removing the product, appropriate aseptic procedures must be observed.

Resterilisation of the product is prohibited.

INSTRUCTIONS FOR USE

Open outer packaging using standard aseptic "no touch" technique to avoid contact with contents. Aseptic personnel only must handle the inner pack.

Do not add any additional substances to the paste. Use only the mixing solutions provided. Using alternative mixing solutions and/or adding other substances to the mixture may alter the setting time significantly. Some substances such as bone marrow and blood will prevent the paste from setting.

MIXING INSTRUCTIONS

BEAD FORMATION

To obtain paste suitable for bead formation, the following steps should be used:-

1. Empty powder into the sterile mixing bowl provided.
2. Add mixing solution from tube provided.
3. Mix thoroughly until a smooth paste is formed (approximately 30 seconds). **DO NOT** over mix.
STIMULAN® Rapid Cure has a working time of 1-2 minutes and will set approximately 4 minutes after mixing.
4. Select the size of bead required and apply a uniform layer of paste onto the bead mat provided. Use the paste applicator to ensure complete filling of each bead cavity.
Note: 10cc STIMULAN® Rapid Cure will fill all of the 3mm or the 4.8mm or the 6mm cavities in the bead mat.
5. Allow paste to set undisturbed for at least 8 minutes after mixing. Flex bead mat to release beads.
6. Before implanting the beads, it is recommended that the structural integrity is checked by compressing the bead between thumb and index finger.

SURGICAL TECHNIQUE

- Accurate diagnosis of any bone cyst is an essential pre-requisite.
- Expose the bone defect. Raise a window of cortical bone with attached muscle or soft tissue where possible. Alternatively raise periosteum without dissection in the extraperiosteal plane. Decortication is to be preferred by sharp osteotome to raise small slivers of bone attached to the deep surface of periosteum.
- In the case of a cyst, evacuate by curette or preferably a high-speed burr. A dental mirror or arthroscope improves clearance of all unwanted tissue by improving visualisation.
- Lavage the prepared bone graft site with a saline solution until a healthy vascular surface is developed on all sides. Beads should not be crushed, shaved or cut.
- Determine the size of the bone void by filling from a syringe of saline. The bone defect should then be gently packed with STIMULAN® Rapid Cure beads or paste.
- Avoid overfilling the bone void or compressing the defect site. Remove any excess paste from the defect site while the material is still workable prior to setting. Do not disrupt the paste while it is setting.
- When implanting the material as a paste, the operating surgeon must ensure a dry (bloodless) field.
- Cover the opening of the bone void with the preserved cortical bone window and/or healthy periosteum or other soft tissues.
- It is of benefit to create a separate micro environment deep to the periosteum. Clean the wound in layers without tension, using standard closure techniques.
- Discard any unused STIMULAN® Rapid Cure material and mixing tools.

Mixing video can be viewed on www.biocomposites.com

For additional information contact your Biocomposites representative or customer service at the address provided.

(EN) ENGLISH

STIMULAN® for use in the presence of Infection

PRODUCT DESCRIPTION

STIMULAN® Rapid Cure is calcium sulfate powder and aqueous mixing solution supplied in pre-measured quantities. When mixed, the resultant paste may be implanted into open voids/gaps of the musculoskeletal system to set in situ or placed in the bead mould mat provided to form beads for implantation. STIMULAN® Rapid Cure is biodegradable and is absorbed by the body in approximately 30 to 60 days. Accessories are provided for the preparation and use of STIMULAN® Rapid Cure. STIMULAN® Rapid Cure, accessories and packaging are not made from natural rubber latex. STIMULAN® is non-metallic and non-magnetic and poses no known hazards in all Magnetic Resonance environments. The characteristics of STIMULAN® Rapid Cure make it suitable as a carrier material for substances such as antibiotics.

INTENDED USES (Not for Canada)

Calcium matrix for bone and soft tissue implantation

INTENDED USES (Canada only)

STIMULAN® Rapid Cure is intended for use as a bone void filler and when mixed with antibiotics, is intended to assist in wound healing and treatment of bacterial infection of the soft tissue surrounding the bone void.

INDICATIONS FOR USE (Not for Canada)

STIMULAN® Rapid Cure is intended for use as a bone / soft tissue void filler for voids or gaps that are not intrinsic to the stability of the bony structure. These osseous defects may be surgically created or created from traumatic injury to the bone.

INDICATIONS FOR USE (Canada Only)

STIMULAN® Rapid Cure is intended for use as a bone void filler for voids or gaps that are not intrinsic to the stability of the bony structure. These osseous defects may be surgically created or created from traumatic injury to the bone.

STIMULAN® Rapid Cure is intended to be gently packed into bony voids or gaps of the skeletal system (i.e. extremities, pelvis, and posterolateral spine). STIMULAN® Rapid Cure provides a bone graft substitute that absorbs and is replaced with bone during the healing process. STIMULAN® Rapid Cure is biodegradable and biocompatible and may be used at an infected site.

To fill voids, defects or gaps created by:

- surgery
- a cyst
- a tumour
- osteomyelitis
- traumatic injury

CONTRAINDICATIONS

- Filling of defects that are intrinsic to the stability of the bony structure
- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative bone disease
- Pregnancy
- Uncooperative patients who can't or won't follow post-operative instructions, including those that abuse drugs or alcohol
- Hypercalcemia

WARNINGS

- Device performance characteristics may be affected by failing to follow indications, contraindications, warnings, precautions or instructions for use
- Concurrent use of other substances may affect bone formation, absorption rate and/or setting time
- Rapid absorption may result in elevated serum calcium levels
- Implantation into patients with a calcium metabolism disorder (including renal compromised patients) may exacerbate symptoms associated with hypercalcemia and/or renal injury and require intervention
- The device has not been studied in paediatric populations
- Do not implant in the cranium where there is risk of contact with cerebrospinal fluid (CSF) or dura mater
- Single use. Do not resterilise. Reuse may result in contamination and/or transmission of disease and/or inconsistent setting characteristics
- Sterile. Do not use beyond expiry date
- Avoid contact with chemicals (e.g. alcohol) in the clinical setting
- Do not implant in patients with known sensitivity or allergic reaction to the device materials

PRECAUTIONS

- This prescription device is for use by suitably qualified Healthcare Professionals only according to its intended use
- It is at the treating physician's discretion to determine the best course of care based on individual patient's needs, taking into account the indications, contraindications, warnings, precautions and possible adverse effects.
- Ensure surrounding bone is healthy and vascular but not actively bleeding at time of insertion
- Do not irrigate the site following implantation
- Beads should not be crushed, shaved or cut to prevent debris restricting pores of the device
- Avoid overfilling or pressurising the defect site (e.g. when injecting) due to risk of migration and/or embolization of fat into the bloodstream
- Ensure adequate tissue coverage and/or containment is achieved as incomplete or inadequate soft tissue coverage may result in migration and/or serous exudate
- Dispose of unused, soiled or explanted product as contaminated clinical waste

POSSIBLE ADVERSE EFFECTS

- Inadequate or incomplete bone formation may occur
- Material sensitivity/allergic reactions may occur
- Heterotopic ossification may occur
- Device migration may occur and result in wound drainage
- Embolism may occur if injecting the device
- Renal failure may occur due to concurrent therapies and/or the patient's medical condition, including implantation into patients who have, or are at risk of, impaired renal function
- Hypercalcemia, infection and seroma are associated risks
- Known risks and potential adverse effects associated with general surgical procedures may also occur

PATIENT INFORMATION

- The patient is to be cautioned to protect the surgical site from unreasonable stresses and to follow the instructions for follow-up care to reduce risk of device migration during the healing period
- The patient is to be made aware of risks and possible adverse effects and to report any such effects, fever or unexpected incidents

STORAGE

Store in a cool, dry place away from direct sunlight and other sources of heat

STERILISATION

- Sterile unless packaging is damaged or opened, check before use
- Indicator dot on packaging turns from yellow to red to indicate product has been exposed to radiation for sterilisation.
Only use product if red
- Double packed to allow aseptic presentation

PREPARATION AND USE

Select size required according to defect or clinical need. Double packed for aseptic presentation. Circulating (non-scrubbed) practitioner to open outer pack and allow the scrubbed practitioner to take the inner pack by grasping or using forceps. Scrubbed practitioner to open inner pack within the sterile field. Prepare STIMULAN® Rapid Cure per surgical requirements. The concurrent use of antibiotics is at the discretion of the physician. Apply STIMULAN® Rapid Cure according to preferred surgical technique.

BEAD FORMATION

1. Empty powder into the sterile mixing bowl provided.
2. Add mixing solution from tube provided. If mixing with antibiotic in liquid form add 6ml per 10cc dose. Refer to warning concurrent use
3. Mix thoroughly until a smooth paste is formed (approximately 30 seconds). DO NOT over mix. STIMULAN® Rapid Cure has a working time of 1-2 minutes and will set approximately 4 minutes after mixing.
4. Select the size of bead required and apply a uniform layer of paste onto the bead mat provided. Use the paste applicator to ensure complete filling of each bead cavity. Note: 10cc STIMULAN® Rapid Cure will fill all of the 3mm or the 4.8mm or the 6mm cavities in the bead mat.
5. Allow paste to set undisturbed for at least 8 minutes after mixing. Flex bead mat to release beads.
6. Before implanting the beads, it is recommended that the structural integrity is checked by compressing the bead between thumb and index finger to confirm the bead has set before implantation.

WARNING

Concurrent use of locally administered antibiotics may affect the setting time, absorption characteristics and / or bone formation. Antibiotics may affect the setting time as exemplified below for 10cc of STIMULAN® Rapid Cure. Only the antibiotics shown can be added and the combination with other antibiotics shall be avoided.

It is the surgeon/healthcare professional's responsibility to give due consideration to the details contained in the medicinal product marketing authorisation in deciding whether it is appropriate for the patient under his/her care. The relevant Summary of Product Characteristics (SmPc) must be consulted. The type and dose of medicinal substance should also be assessed according to the individual patient's clinical circumstances.

Data based on 10cc Stimulan Rapid Cure for a given quantity of substance. Supplier shown in parenthesis, other sources may vary.

For additional information contact your Biocomposites representative or distributor.

Single powders added to Stimulan powder plus mixing solution provided

Substance (supplier)	Quantity	Expected Set Time (min)
Vancomycin (Flynn Pharma)	1000mg	2-5 min

Single Liquids added to Stimulan powder instead of mixing solution provided

Substance (supplier)	Quantity	Expected Set Time (min)
Gentamicin (Amdipharm)	240mg in 6ml solution	3-5 min
Tobramycin (Alloga)	240mg in 6ml solution	5-15 min

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Biocomposites Ltd, Keele Science Park, Keele, Staffordshire, ST5 5NL, England
Biocomposites Inc, PO Box 2692, Wilmington, NC 28402, USA