I.CO.S.[®]

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NEWDEAL[®] as the manufacturer of this device, does not practice medicine and does not recommend this or any other surgical technique for use on a specific patient.

The surgeon who performs any implant procedure is responsible for determining and using the appropriate techniques for implanting the device in each patient.

Stage 1 POSITIONING OF THE PIN:

A 1.6 mm diameter K-wire, 150 mm long (115 116), is introduced into the site where the I.CO.S screw is to be implanted. The K-wire is implanted by means of a drilling guide (119 506) inserted in the protecting sleeve (119 501). Failure to use the drilling guide and protecting sleeve when inserting the K-wire may result in soft tissues entrapment. Care should be taken when K-wire is inserted to respect neighbouring blood vessels, nerves and tissues. The depth of insertion of the K-wire depends upon the surgeon's will to achieve a mono or bi-cortical osteosynthesis.

Stage 2 DETERMINATION OF THE LENGTH OF SCREW TO BE IMPLANTED:

At insertion of the k-wire, the appropriate length of screw is determined with a measurer (119 510). The measurer is placed so that the thin side is aligned with the base of the K-wire. Care should be taken to hold the measurer in perfect alignment with the visible part of the K-wire. The appropriate screw length is then read directly from the measurer (graduated from 0 to 60 mm), on the side marked «K-wire 150 mm».


The surgeon may decide to add or deduct a few millimeters to / from the length read from the ruler. For example, a few millimeters may be deducted to ensure that the second cortex will not be penetrated.

IMPORTANT:

The measurement taken from the ruler / K-wire is directly related to the positioning of the K-wire and cannot therefore be considered ideal unless the K-wire is correctly positioned. The positioning of the K-wire can be validated by means of peroperative X-ray or fluoroscopy.

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The window on the back of the ruler (119 510) can be used to confirm the length of the screw to be implanted. The head of the I.CO.S.[®] screw is positioned on the window; the distal part of the screw directed towards the thin part of the measurer. The mark 0 on the internal part of the measurer also allows for this control.

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Stage 3

USE OF DRILLS:

Once the 1.6 mm K-wire is in place, the drills are used to prepare the insertion of the screw.

3.1 • The short cannulated drill (119 522) attached on the quick coupling (129 710) is used to prepare the location of the screw head. The protecting sleeve (119 501) must be used during this stage.

The short drill will be inserted in the sleeve before drilling. This allows to limit the depth of reaming thanks to marks on the drill.

When the drill is used perpendicular to the bone, the mark 90° must be flush to the extremity of the sleeve.

When the drill has a 45° angulation, the mark 45° will be aligned with the extremity of the sleeve. It is only when there is a 30° degree orientation that the drill will be fully engaged in the protection sleeve.

This allows for reaming only to a depth of 5 mm sufficient to prepare the location of the screw head.

3.2 • For all the screw sizes – unless a mono-cortical fixation is required, the cannulated drill (119 516) is used to finalise the preparation stage. This drill is introduced over the pin and drilled to a depth that corresponds to the length of the screw to be implanted.

CAUTION : Drills should never be used without the matching pin.

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NEWDEAL[®] as the manufacturer of this device, does not practice medicine and does not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and using the appropriate techniques for implanting the device in each patient.

Stage 1 **POSITIONING OF THE PIN:**

A 2.5 mm diameter graduated K-wire, 200 mm long (115 225) is introduced into the site where the I.CO.S. screw diam 6.5 mm is to be implanted. The K-wire is positioned by means of drilling guides and soft tissues protector (119541, 119551 and 119552). Failure to use the drilling guide and protecting sleeve when inserting the K-wire may result in soft tissues entrapment. Care should be taken when K-wire is inserted to respect neighbouring blood vessels, nerves and tissues. The depth of insertion of the k-wire depends upon the surgeon's will to achieve a mono or bi-cortical osteosynthesis.

Stage 2 **DETERMINATION OF THE LENGTH OF SCREW TO BE IMPLANTED:**

At insertion of the k-wire, the appropriate length of screw is determined from a measurer (119 510). The measurer is placed so that the thin side is aligned with the base of the k-wire. Care should be taken to hold the measurer in perfect alignment with the visible part of the k-wire. The appropriate screw length is then read directly from the measurer (graduated from 0 to 90 mm), on the side marked «K-wire 200 mm».


The surgeon may decide to add or deduct a few millimeters to/from the length read from the ruler. For example, a few millimeters may be deducted to ensure that the second cortex will not be penetrated.

IMPORTANT:

The measurement taken from the ruler/k-wire is directly related to the positioning of the k-wire and cannot therefore be considered ideal unless the K-wire is correctly positioned. The positioning of the K-wire can be validated by means of peroperative X-ray or fluoroscopy.

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The superior part of the measurer (119 510) allows for the control of the screw to be implanted. The head of the I.CO.S.* screw is positioned on the top of it, the distal part directed towards the thin part of the measurer.

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Stage 3

USE OF DRILLS:

Once the 2.5 mm diameter K-wire is in place, the drills are used to prepare the insertion of the screw.

3.1 • The short cannulated drill (119 546) attached on the quick coupling (129 710) is used to prepare the location of the screw head. The external protecting sleeve (119 552) must be used during this stage.

The short drill will be inserted in the sleeve before drilling. This allows to limit the depth of reaming thanks to marks on the drill.

When the drill is used perpendicular to the bone, the mark 90° must be flush to the extremity of the sleeve.

When the drill has a 45° angulation, the mark 45° will be aligned with the extremity of the sleeve. It is only when there is a 30° degree orientation that the drill will be fully engaged in the protection sleeve.


This allows for reaming only to a depth of 10 mm sufficient to prepare the location of the screw head.

3.2 • The countersinking for the outer thread thus prepared, the long drill (119 571) is used to prepare the distal cortex for introduction of the screw. The use of this long drill is not recommended unless a bi-cortical fixation is required. This long drill is introduced over the K-wire left in place.

It is graduated so that the drilled depth can be read directly. The visible number the nearest to the point of entry of the drill into the bone is the length of the buried part of the drill.

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Stage 4

USE OF THE TAP:

With the K-wire and the protecting sleeve still in place, the tap (119 555 / 119 525) is now introduced to prepare the proximal cortex. The tap is then screwed into the countersink area intended for the screw head, then unscrewed, taking care not to damage the prepared threads. The in-depth of the tap should correspond to the angulation confirmed during the previous step.

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Stage 5

SCREW INSERTION:

With the K-wire still in place, the I.CO.S. screw is selected and introduced by means of a cannulated screwdriver tip (119 562 / 119 532) attached on the quick coupling, or manually with the internal screwdriver (119 561 / 119 531). The extremity of the head of the screw should be flush to the bone cortex. At this point static compression has been achieved.

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NB: Care should be taken that the external part is completely locked before insertion on the K-wire.

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Stage 6

OPTIONAL DYNAMIC COMPRESSION:

6.1 • The static compression effect can be increased by a continuous screwing movement of the entire screw. However, in order to maintain effective compression and adequate stability the head of the screw should remain in the cortical bone.

6.2 • Additional (optional) dynamic compression is achieved by advancing the outer thread of the screw independently of its main body. For this, the external screwdriver (119 566 / 119 536) is used, allowing for screwing the external threading of the screw head. The body (internal part) of the screw will not be moved.

Stage 7

OPTIONAL TRANSLATION OF THE SCREW BODY:

Once dynamic compression has been used, this may leave the proximal end of the screw out of the cortex. By using the internal screwdriver (119 561 / 119 531) the screw can be further implanted. This does not in any way reduce the effect of the compression because the distal and proximal threads of the body of the screw are identical.

Only after this stage can the guide wire be removed.

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

Mèche épaulée
Short cannulated drill

Mèche longue
Long cannulated drill

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INSTRUCTIONS FOR USE

NON STERILE IMPLANTS FOR FOOT SURGERY • SINGLE USE

In accordance with EEC directive 93/42 relative to medical devices, this product must be handled and/or implanted by WELL-TRAINED, QUALIFIED PERSONS, AWARE OF THESE DIRECTIONS FOR USE.

1 - Description of the medical devices:

The implants - delivered non sterile - are:

- Compression screws existing in different diameters and lengths;
- They are made out of Titanium alloy within the frame of the standard NF ISO 5832-3 and ASTM F136.

2 - Indications:

The I.C.O.S. Screw is indicated for fixation of bone fractures or for bone reconstruction.

Examples include:

- Fixation of bone fragments, in long bones or small bones fractures;
- Fractures management in the foot or hand;
- Arthrodesis in hand, foot or ankle surgery;
- Mono or Bi-cortical osteotomies in the foot or hand or in long bones;
- Treatment of inferior tibio fibular diastasis.

The size of the chosen screw should be adapted to the specific indication.

3 - Contraindications:

The implant should not be used in a patient who has currently, or who has a history of:

- Local or systemic acute or chronic inflammation;
- Active infection or inflammation;
- Suspected or documented metal allergy or intolerance.

4 - Warnings:

Serious post-operative complications may occur from use of the implant in a patient who:

- Lacks good general physical condition;
- Has severe osteoporosis;
- Demonstrates physiologic or anatomic anomalies;
- Has immunological responses, sensitization, or hypersensitivity to foreign materials;
- Systemic or metabolic disorders.

5 - Precautions for use:

Physician must determine if implant is appropriate for patients who have any of the following conditions:

- Drug and/or alcohol and/or smoke addiction and/or abuse;
- Infectious disease;
- Malignancy;
- Local bone tumors;
- Systemic or metabolic disorders or replacement;
- Compromised wound healing;
- Obesity;
- Demonstrated psychological instability, displayed a lack of understanding, inappropriate motivation, or attitude;
- Unwillingness to accept the possibility of multiple surgeries for revision or replacement;
- Lacks an understanding that a metallic implant is not as strong as normal healthy bone and will bend, loosen, or fracture if excessive demand is placed on it;
- Lacks an understanding that their preoperative capacity may not be fully recovered even after successful implantation.

Knowledge of surgical techniques, proper reduction, selection and placement of implants, and post-operative patient management are considerations essential to a successful outcome.

Criteria for patient selection is the responsibility of the surgeon. Information contained within this document should be taken into consideration during the selection process. Recognition of the appropriate indications and contraindications and the selection of the proper surgical procedures and techniques determined to be best for the patient are the responsibility of the surgeon. Each surgeon must evaluate the appropriateness of the procedure and instruments used during the procedure based on his or her own training and experience.

The surgeon should discuss with the patient prior to surgery possible risks, precautions, warnings, consequences, complications, and adverse reactions associated with the surgical procedure and implantation of the device.

Each patient must be evaluated by the surgeon to determine the specific risk/benefit relationship in light of the patient's condition and the surgeon's practice, training, experience, and knowledge of the related medical literature.

Complications with the use of compression screws have been reported in the medical literature. Any patient undergoing a surgical procedure is subject to intra-operative and post-operative complications. Each patient's tolerance to surgery, medication, and implantation of a foreign object may be different.

Possible risks, adverse reactions, and complications associated with surgery and the use of the compression screws should be discussed with and understood by the patient prior to surgery. The implant is composed of titanium alloy materials; therefore, it is subject to possible reactions and complications, including those listed herein. The patient should not be led to unrealistic.

expectations as to the performance or results that the surgery and implant can provide. The patient should be informed that the life expectancy of the device is unpredictable once implanted, and that successful results cannot be guaranteed.

IT IS THE RESPONSIBILITY OF THE SURGEON TO PROVIDE THE PATIENT WITH INFORMATION PRIOR TO SURGERY.

Complications may include but are not limited to:

- Pain, discomfort, or abnormal sensations due to presence of the implant;
- Bending, loosening, and/or breakage, which could make removal impracticable or difficult;
- Risk of additional injury from post-operative trauma;
- Migration of the implant position or implant material resulting in injury;
- Bone loss due to stress shielding;

Side effects may include but are not limited to:

- Infections;
- Hematoma;
- Allergy;
- Thrombosis;
- Bone non union or delayed union.

Adverse effects may necessitate re-operation, revision or removal surgery, arthrodesis of the involved joint, and /or amputation of the limb.

Implant removal should be followed by adequate postoperative management to avoid fracture or re-fracture.

Interference risks during medical imaging:

MRI/SCANNER: ask the patient to systemically mention that he/she was implanted with a metallic device.

6 - Instrucons for reprocessing:

This product is sold non sterile.

Check the integrity of the packaging and labeling before opening the packing.

Remove all the products from their packaging prior to sterilization.

All products should be cleaned, decontaminated, and sterilized before use.

Always immediately clean and decontaminate all devices that have been soiled.

Repeated reprocessing has little effect on these products.

Preparation: Double instruments (ex. Internal screwdriver and associated external screwdriver) should be separated prior to cleaning.

Cleaning: Cleaning can be performed manually, automatically or ultrasonically in accordance with the specifications designated by the manufacturer of the hospital's equipment.

Manual cleaning: Manual cleaning consists of using aldehyde free cleaners (neutral or alkaline), applied with a soft brush, taking special care to threaded parts and parts difficult to reach.

Note: Certain solutions such as those containing bleach or formalin may damage the devices, and they must not be used.

Use of metallic brushes or other abrasive products is also forbidden.

Cleaning should be immediately followed by profusely rinsing with deionized water. Check that water flows out the cannulated parts.

Automatic cleaning: Automatic cleaning is performed in a cleaning/disinfecting machine using neutral cleaners, with a cleaning cycle of 5 minutes minimum and a rinsing cycle of 3 minutes.

Check the complete removal of visible dirt, especially in the cannulated parts.

If necessary, repeat the full process or proceed to a manual cleaning.

Disinfection: If an automatic cleaning is used, final rinsing at 95°C during 10 minutes can be performed.

Drying: Drying temperature should not exceed 95°C.

Controls, servicing and tests: No specific requirements.

The implants are single use. They should therefore never be re-used.

Packaging: No specific requirements.

Sterilisation:

Newdeal's implants and instruments are recommended to be sterilized by the steam autoclaving procedure regularly used in the hospital.

The following two methods have been validated by the manufacturer and can thus be used:

Method: steam	Method: steam
Cycle: wrapped gravity	Cycle: wrapped gravity
Temperature: 132°C	Temperature: 134°C
Exposure time: 45 minutes	Exposure time: 18 minutes

Other sterilization method and cycles may also be used. However, individuals or hospitals not using the recommended method are advised to validate the alternative method using appropriate laboratory techniques. EtO sterilization or cold sterilization techniques are not recommended.

7 - Use of the implant:

The surgeon must use the instrumentations recommended in accordance with the operative technique available from the manufacturer. The medical device must be used in compliance with the use of the profession and the standard of art. Do not attempt a surgical procedure with faulty, damaged or suspect instruments or implants. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.

Opening of the instruments set must be done according to aseptic condition.

When handling the implants, avoid any contact with other material or tools which may damage the implant surface.

Under no circumstances should the implant be modified.

8 - Re-use of the implants:

Orthopaedic implants already implanted must never be re-used. The company accepts no responsibility for such re-use.

9 - Re-sterilization of non implanted products:

Re-sterilization is only allowed for non implanted products. Such non implanted products can be sterilized several times in the same conditions as those described above.

10 - Preventative actions for the patient to avoid post-operative complications:

- Avoid extreme position such as flexion-extension;
- wear orthopaedic shoes according to the surgeon's prescription;
- receive prompt medical attention for any infection that could occur, whether at the operated-member level or elsewhere in the body.

11 - Storage: Store in dry place.

12 - Product information disclosure / Liability:

Newdeal, an Integra LifeSciences Company, has exercised reasonable care in the selection of materials and the manufacture of these products. Newdeal excludes all warranties, whether expressed or implied, including but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Newdeal shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of this product. Newdeal neither assumes nor authorizes any person to assume for it any other or additional liability or responsibility in connection with these products. Newdeal intends that this device should be used only by physicians having received proper training in orthopedic surgery technique for use of the device.

WARNING: Federal law (USA) restricts this device to sale by or on the order of a physician.

WARNING: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

INFORMATION: Should any information regarding the products or their uses be required, please contact your representative or distributor or directly contact the manufacturer.

I.CO.S. • DIAM. 4.0 MM

REFERENCE	DESCRIPTION
105 426	L. 26 mm
105 428	L. 28 mm
105 430	L. 30 mm
105 432	L. 32 mm
105 434	L. 34 mm
105 436	L. 36 mm
105 438	L. 38 mm
105 440	L. 40 mm
105 445	L. 45 mm
105 450	L. 50 mm
105 455	L. 55 mm
105 460	L. 60 mm

I.CO.S. • DIAM. 6.5 MM.

REFERENCE	DESCRIPTION
105 640	L. 40 mm
105 645	L. 45 mm
105 650	L. 50 mm
105 655	L. 55 mm
105 660	L. 60 mm
105 665	L. 65 mm
105 670	L. 70 mm
105 675	L. 75 mm
105 680	L. 80 mm
105 685	L. 85 mm
105 690	L. 90 mm

INSTRUMENTS ASSOCIES • ASSOCIATED INSTRUMENTS

REFERENCE	DESCRIPTION
119 501	CANON DE PROTECTION PROTECTION SLEEVE
119 506	GUIDE DRILL GUIDE
119 510	MESUREUR MEASURER
119 516	MÈCHE DROITE LONG CANNULATED DRILL
119 522	MÈCHE ÉPAULÉE SHORT CANNULATED DRILL
129 710	SYSTÈME DE FIXATION RAPIDE QUICK-COUPLING
119 525	TARAUD TAP
119 531	TOURNEVIS INTERNE INTERNAL SCREWDRIVER
119 532	EMBOUT DE TOURNEVIS SCREWDRIVER TIP
119 536	TOURNEVIS EXTERNE EXTERNAL SCREWDRIVER
119 580	TREFINE TREPINE

INSTRUMENTS ASSOCIES • ASSOCIATED INSTRUMENTS

REFERENCE	DESCRIPTION
119 510	MESUREUR MEASURER
119 541	GUIDE DRILL GUIDE
119 551	CANON DE PROTECTION INTERNE INTERNAL PROTECTION SLEEVE
119 552	CANON DE PROTECTION EXTERNE EXTERNAL PROTECTION SLEEVE
119 571	MÈCHE LONGUE LONG CANNULATED DRILL
119 546	MÈCHE ÉPAULÉE SHORT CANNULATED DRILL
129 710	SYSTÈME DE FIXATION RAPIDE QUICK-COUPLING
119 555	TARAUD TAP
119 561	TOURNEVIS INTERNE INTERNAL SCREWDRIVER
119 562	EMBOUT DE TOURNEVIS SCREWDRIVER TIP
119 566	TOURNEVIS EXTERNE EXTERNAL SCREWDRIVER
119 590	TREFINE TREPINE

REFERENCE	DESCRIPTION
115 116	BROCHE DE KIRSCHNER K-WIRE DIAM. 1.6 MM / L. 150 MM
115 225	BROCHE DE KIRSCHNER K-WIRE DIAM. 2.5 MM / L. 200 MM
119 974	CONTENEUR DE STÉRILISATION I.CO.S.® ø 4,0 MM STERILIZATION CONTAINER I.CO.S.® DIAM 4.0 MM
119 976	CONTENEUR DE STÉRILISATION I.CO.S.® ø 6,5 MM STERILIZATION CONTAINER I.CO.S.® DIAM 4.5 MM

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- The products are manufactured and referenced within the frame of the standards in force.
- Implantation procedures are described in the surgical technique.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- **WARNING** : Federal law (USA) restricts this device to sale by or on the order of a physician.

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NEWDEAL SAS