

TenaTac[®]

LTT001.00 rev.2019-10

EN	Instructions for use	4
NL	Gebruiksaanwijzing	6
DA	Brugsvejledning	8
DE	Gebrauchsinformation	10
FR	Notice d'utilisation	12
ES	Instrucciones de uso	14
PT	Instrucoes para utilizacco	16
TR	Kullunim talimati	18
PL	Instrukcje użytkowania	20
IT	Spugna di gelatina riassorbibile	22
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	24
CZ	Návod k použití houby	26
SK	Návod na použitie špongie	28
RO	Instrucțiuni pentru utilizare	30
SR	Uputstvo za korišćenje	32
HR	Uputa za uporabu	34

Instructions for use TenaTac®

1. Product description

The TenaTac® gelatine sponge is made of 100% purified porcine gelatine and possesses a unique surface modification for improved adhesion to wounds. It is intended for the use as an absorbable haemostat in various fields of surgery, where the significantly improved adhesion to wounds ensures that the haemostat stays in place. It is insoluble in water and fully absorbable. Its porous structure and interstices enable the product to absorb more than 35 times its own weight of blood and fluids. The absorptive capacity of TenaTac® is a function of its physical size, increasing as the amount of gelatine increases. The unique surface of TenaTac® significantly increases tissue adhesion relative to plain gelatine sponge. The surface modification allows TenaTac® to be more flexible than plain gelatine sponge. The improved adhesion has been achieved without having to incorporate expensive coagulants or reactive chemicals into the sponge. However, because of its neutral character and nature, it may be used as a carrier for any drug (the product literature for the applicable drug should be consulted for complete prescribing information). It is supplied in several sizes. It is packaged individually in (double) blister, pouch or container packaging, sterilized by gamma rays, and is for single use only.

2. Indications

- 2.1 General indications: the TenaTac® gelatine sponge is used during and after surgical procedures to achieve local haemostasis by tamponade effect.
- 2.2 For ENT surgery the use of gelatine sponge has been proven to be effective in supporting and keeping the fascia or perichondrium in the middle ear (dry application) in place or in the outer ear canal after tympanoplasty.
- 2.3 For rectal surgery, haemorrhoid operations and gynaecology TenaTac® gelatine sponge can be used to control bleeding.

3. Contra-indications

- Allergies to porcine products.
- The gelatine sponge should not be used in closures of skin incisions, as it may interfere with

wound healing.

- The product should not be used without antibiotics in infected wounds.

4. Side effects

There have been no reported adverse effects for gelatine sponge when the product is used correctly according to these Instructions for Use and when over-packing is not applied. Formation of tissue granulation during ENT/middle ear procedures has been reported in at least one animal study utilising the base gelatin sponge.

5. Directions for use

TenaTac® is flexible and ready-to-use in dry form; no wetting is required. For special indications, such as the carriage of an active compound, the adhered sponge can be soaked with a saline solution. Moistening the sponge prior to application will increase flexibility but can reduce overall adhesion.

When used dry, the packing is readily opened and the sponge is cut into the desired size, usually a slightly larger than the bleeding area and is lightly compressed into position. The sponges must be applied to a bleeding area under light pressure with a moist open-weave gauze for one or two minutes until the bleeding stops. It is possible to retreat a bleeding site with further applications of TenaTac® if needed.

If bleeding is controlled after (possible multiple) applications of the gelatine sponge, the material can be left in situ. The excess should be removed.

Depending on the method of use and the operated site, the total absorption time of the sponge varies from four to six weeks. When applied to bleeding nasal, rectal or vaginal mucosa, it liquefies within 2 to 5 days. It should be noted that the absorption is a process where the volume as well as the absorption characteristics continuously decrease over time.

If moistening with saline is required, TenaTac® should be soaked in the solution, then withdrawn, squeezed thoroughly between gloved fingers to expel air bubbles present in the interstices, replaced in saline, and kept there until needed.

5.1 ENT applications:

The gelatine sponge is cut as needed to fit the otic, nasal or oral cavity or defect and inserted to support and separate tissues to prevent adhesions and/or to control bleeding by tamponade effect. The excess should be removed.

Important notice!

In cases where the gelatine sponge has been compressed (either by the manufacturer or end-user), the gelatine sponge may expand with the absorption of fluids (the sponge will obtain its original shape). Where possible TenaTac® should be removed after usage and bleeding has stopped in radical cavities, laminectomy procedures, around or in proximity to foramina in bone, areas of bony confine, the spinal cord and/or the optic nerve and chiasm or closed tissue spaces with presence of bone. This might lead to unintended pressure on neighbouring structures which may result in pain for the patient or might create the potential for nerve damage.

The over-packing of TenaTac®, particularly within bony cavities, should be avoided, since recovering to its initial volume may interfere with normal function and/or could cause possible or eventual compression necrosis of surrounding tissue. Formation of tissue granulation has been observed during ENT / middle ear procedures with plain gelatine sponges.

In cases of postoperative infection, re-operation may be necessary to remove infected material and allow drainage.

Warnings

Do not use when the package has been opened or has been damaged.

The absorbable gelatine TenaTac® sponge is a single-use product which is not suitable for re-sterilization.

6. General storage and handling

The product should be stored in its sealed protective outer packaging. Prior to opening, the packaging ensuring the integrity of the sterile barrier has

to be inspected. Open the package just prior to use in a patient.

6.1 Packaging and sterility

The products are packaged sterile.

The packaging consists of:

- Storage and transport packaging (secondary packaging)
- Sterile packaging (primary packaging)

The packaging is subject to and conforms to European & International legislation and other applicable standards. The packaging protects the product from external influences and guarantees its sterility during storage.

6.2 Handling of the sterile packaging

Remove the product from the sterile package using acceptable aseptic technique.

6.3 Sterilization

The products are sterilized by Gamma radiation.

6.4 Storage

The products should be stored in a dry place in its protective packaging at room temperatures. Avoid contact with direct sunlight. The expiry date of TenaTac® is printed on the pack. Do not use TenaTac® after this date.

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