



# KellMed Ltd.

KELLMED SALES LTD IS A SPECIALIST SUPPLIER  
OF STERILE DISPOSABLE INSTRUMENTS FOR  
ORTHOPAEDIC AND TRAUMA SURGERY.

---





- We offer a core product range of K Wires, Drill Bits, Steinmann Pins, Guide Wires and Bone Taps – all sterile packed for single use.
- All our products are CE marked and are packed and sterilised in the UK.
- KellMed works to a strict code of practice to certification ISO 13485 and ISO 9001. We have been supplying NHS and private hospitals across the UK for over ten years. Our focus on efficiency at every stage of the supply chain ensures that our customers receive top quality products at highly competitive prices.



In the late 1990s we became interested in the impact of re-use on the quality of orthopaedic drill bits. Recycling (cleaning, sterilising, re-using) was standard practice, but there was no way to reliably identify or sharpen drill bits that had become worn with use. As a result, blunted bits could be passed to the surgical team, potentially impacting clinical outcomes. In a parallel development, concerns emerged over the potential for HIV and variant CJD transmission by reused surgical devices. Therefore, in 2002 we formed a company to produce and deliver high quality sharp products, sterile packed for single use, seeing this as a practical solution to both problems.

Our strategic focus on high-quality drill bits, taps, pins and wires helps us maintain competitive prices and the inventory required to meet your requirement for rapid response and delivery. We are also able to produce small quantities of specialised products to meet specific needs.

# drill bits

DRILL BITS WITH QUICK FIT CONNECTION

DRILL BITS WITH PLAIN SHANK



## All our drills are:

CE marked and manufactured to international standards

Individually packed in tough Tyvek® pouches

Sterile with a shelf life of five years

Supplied with three traceability labels on outer pouch of each item

Available in other sizes and fittings – please let us know your requirements

drill bits

drill bits with quick fit connection



Drill bits with AO-type connection

Description	Product Code
1.1 mm x 60 mm drill bit, quick fit	912.060s
1.5 mm x 85 mm drill bit, quick fit	913.085s
1.8 mm x 100 mm drill bit, quick fit	920.100s
2.0 mm x 100 mm drill bit, quick fit	914.100s
2.0 mm x 125 mm drill bit, quick fit	914.125s
2.4 mm x 100 mm drill bit, quick fit	921.100s
2.5 mm x 110 mm drill bit, quick fit	915.110s
2.5 mm x 125 mm drill bit, quick fit	915.125s
2.7 mm x 100 mm drill bit, quick fit	916.100s
2.7 mm x 125 mm drill bit, quick fit	916.125s
3.2 mm x 125 mm drill bit, quick fit	917.125s
3.2 mm x 145 mm drill bit, quick fit	917.145s
3.2 mm x 195 mm drill bit, quick fit	917.195s
3.5 mm x 110 mm drill bit, quick fit	918.110s
3.5 mm x 145 mm drill bit, quick fit	918.145s
3.5 mm x 195 mm drill bit, quick fit	918.195s
4.5 mm x 145 mm drill bit, quick fit	919.145s
4.5 mm x 195 mm drill bit, quick fit	919.195s

**All our drills are:**

- CE marked and manufactured to international standards
- Individually packed in tough Tyvek® pouches
- Sterile with a shelf life of five years
- Supplied with three traceability labels on outer pouch of each item
- Available in other sizes and fittings - please let us know your requirements

## drill bits with plain shank



## Drill bits with plain shank

Description	Product Code
1.1 mm x 45 mm drill bit, plain shank	930.045s
2.0 mm x 85 mm drill bit, plain shank	932.085s
1.5 mm x 45 mm drill bit, plain shank	931.070s
2.5 mm x 95 mm drill bit, plain shank	933.095s
2.7 mm x 100 mm drill bit, plain shank	934.100s
2.7 mm x 85 mm drill bit, plain shank	934.085s
3.2 mm x 106 mm drill bit, plain shank	935.106s
3.2 mm x 130 mm drill bit, plain shank	935.130s
3.5 mm x 180 mm drill bit, plain shank	936.180s
4.0 mm x 120 mm drill bit, plain shank	939.120s
3.5 mm x 112 mm drill bit, plain shank	936.112s
4.0 mm x 180 mm drill bit, plain shank	939.180s
4.0 mm x 120 mm drill bit, plain shank	939.120s
4.5 mm x 180 mm drill bit, plain shank	937.180s
6.5 mm x 150 mm drill bit, plain shank	938.150s


**All our drills are:**

- CE marked and manufactured to international standards
- Individually packed in tough Tyvek® pouches
- Sterile with a shelf life of five years
- Supplied with three traceability labels on outer pouch of each item
- Available in other sizes and fittings - please let us know your requirements

**KellMed**

# wires & pins

K WIRES, SINGLE TROCAR  
K WIRES, DOUBLE TROCAR  
CUSTOM MADE K WIRES  
STEINMANN PINS  
GUIDE WIRES



## **All our wires are:**

CE marked and manufactured to international standards

Individually packed in tough Tyvek® pouches

Sterile with a shelf life of five years

Supplied with three traceability labels on outer pouch of each item

Available in other sizes and fittings – please let us know your requirements



The single trocar K wires on this list are the most popular sizes in our range. They are all held in stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know, we may well have it in stock. If not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.

wires & pins



## k wires, single trocar

Description	Product Code
<b>0.8 mm diameter</b>	
0.8 mm x 80 mm K wire, single trocar	319.080s
0.8 mm x 150 mm K wire, single trocar	319.150s
<b>0.9 mm diameter</b>	
0.9 mm x 150 mm K wire, single trocar	362.150s
0.9 mm x 125 mm K wire, single trocar	362.125s
<b>1.0 mm diameter</b>	
1.0 mm x 100 mm K wire, single trocar	320.100s
1.0 mm x 150 mm K wire, single trocar	320.150s
1.0 mm x 200 mm K wire, single trocar	320.200s
<b>1.1 mm diameter</b>	
1.1 mm x 150 mm K wire, single trocar	363.150s
1.1 mm x 200 mm K wire, single trocar	363.200s
<b>1.2 mm diameter</b>	
1.2 mm x 150 mm K wire, single trocar	321.150s
1.2 mm x 305 mm K wire, single trocar	321.305s
<b>1.25 mm diameter</b>	
1.25 mm x 150 mm K wire, single trocar	346.150s
1.25 mm x 150 mm K wire, single trocar, threaded tip	347.150s
<b>1.4 mm diameter</b>	
1.4 mm x 150 mm K wire, single trocar	322.150s
<b>1.5 mm diameter</b>	
1.5 mm x 150 mm K wire, single trocar	364.150s
1.5 mm x 200 mm K wire, single trocar	364.200s
<b>1.6 mm diameter</b>	
1.6 mm x 150 mm K wire, single trocar	323.150s
1.6 mm x 150 mm K wire, single trocar, threaded tip	348.150s
1.6 mm x 200 mm K wire, single trocar	323.200s
<b>1.8 mm diameter</b>	
1.8m x 150 mm K wire, single trocar	325.150s
<b>2 mm diameter</b>	
2.0 mm x 75 mm K wire, single trocar	324.075s
2.0 mm x 150 mm K wire, single trocar	324.150s
2.0 mm x 200 mm K wire, single trocar	324.200s
<b>2.5 mm diameter</b>	
2.5 mm x 230 mm K wire, single trocar, threaded tip	353.230s
2.5 mm x 260 mm K wire, single trocar, threaded tip	353.260s



The double trocar K wires on this list are the most popular sizes in our range. They are all held in stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know, we may have it in stock. If not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.

wires & pins

## k wires, double trocar

Description	Product Code
<b>0.8 mm diameter</b>	
0.8 mm x 150 mm K wire, double trocar	338.150s
<b>0.9 mm diameter</b>	
0.9 mm x 75 mm K wire, double trocar	370.075s
0.9 mm x 100 mm K wire, double trocar	370.100s
0.9 mm x 125 mm K wire, double trocar	370.125s
0.9 mm x 150 mm K wire, double trocar	370.150s
<b>1.0 mm diameter</b>	
1.0 mm x 125 mm K wire, double trocar	339.125s
1.0 mm x 150 mm K wire, double trocar	339.150s
<b>1.1 mm diameter</b>	
1.1 mm x 75 mm, K wire, double trocar	371.075s
1.1 mm x 100 mm, K wire, double trocar	371.100s
1.1 mm x 125 mm, K wire, double trocar	371.125s
1.1 mm x 150 mm, K wire, double trocar	371.150s
<b>1.2 mm diameter</b>	
1.2 mm x 100 mm K wire, double trocar	340.100s
1.2 mm x 150 mm, K wire, double trocar	340.150s
1.2 mm x 200 mm K wire, double trocar	340.200s
1.2 mm x 230 mm K wire, double trocar	340.230s
<b>1.25 mm diameter</b>	
1.25 mm x 150 mm, K wire, double trocar	345.150s
<b>1.4 mm diameter</b>	
1.4 mm x 100 mm, K wire, double trocar	341.100s
1.4 mm x 150 mm, K wire, double trocar	341.150s
<b>1.5 mm diameter</b>	
1.5 mm x 150 mm, K wire, double trocar	372.150s
1.5 mm x 300 mm, K wire, double trocar	372.300s
<b>1.6 mm diameter</b>	
1.6 mm x 75 mm, K wire, double trocar	342.075s
1.6 mm x 100 mm, K wire, double trocar	342.100s
1.6 mm x 125 mm, K wire, double trocar	342.125s
1.6 mm x 150 mm, K wire, double trocar	342.150s
1.6 mm x 200 mm, K wire, double trocar	342.200s
<b>1.8 mm diameter</b>	
1.8 mm x 150 mm, K wire, double trocar	343.150s
<b>2.0 mm diameter</b>	
2.0 mm x 125 mm K wire, double trocar	344.125s
2.0 mm x 150 mm, K wire, double trocar	344.150s
2.0 mm x 200 mm K wire, double trocar	344.200s

wires & pins

custom made k wires



Description	Single Trocar	Double Trocar
0.9 mm K wire	362.XXXs	370.XXXs
1.0 mm K wire	320.XXXs	339.XXXs
1.1 mm K wire	363.XXXs	371.XXXs
1.2 mm K wire	321.XXXs	340.XXXs
1.25 mm K wire	346.XXXs	345.XXXs
1.4 mm K wire	322.XXXs	341.XXXs
1.5 mm K wire	364.XXXs	372.XXXs
1.6 mm K wire	323.XXXs	342.XXXs
1.8 mm K wire	325.XXXs	343.XXXs
2.0 mm K wire	324.XXXs	344.XXXs
2.5 mm K wire	335.XXXs	337.XXXs

**Order wires to any length:**

- Add the length (mm) to the main code given in the grid above.  
*For example, all types of 0.9 mm K wire, single trocar, sterile have the code 362.XXXs.*

**Therefore:**

- 0.9 mm K wire, single trocar, length 50 mm, sterile..... 362.050s
- 0.9 mm K wire, single trocar, length 150 mm, sterile..... 362.150s

**All custom made sterile packed wires are:**

- Individually packed in tough Tyvek® backed polyethylene pouches.
- Sterile with a shelf life of five years
- Supplied with three traceability labels per wire
- Available in 6-8 weeks from initial order, after which we will maintain stock ready for rapid delivery.

wires & pins

steinmann pins



Description

Product Code

**Steinmann pins, trocar/triangle**

3.5 mm x 200 mm Steinmann pin, trocar/ triangle	313.200s
4.0 mm x 200 mm Steinmann pin, trocar/ triangle	306.200s
4.0 mm x 175 mm Steinmann pin, trocar/ triangle	306.175s
4.0 mm x 225 mm Steinmann pin, trocar/ triangle	306.225s
4.0 mm x 254 mm Steinmann pin, trocar/ triangle	306.254s
4.8 mm x 175 mm Steinmann pin, trocar/ triangle	305.175s

**Steinmann pins, trocar/triangle, fully threaded**

2.4 mm x 230 mm Steinmann Pin, trocar/ triangle, fully threaded	317.230s
2.8 mm x 230 mm Steinmann Pin, trocar/ triangle, fully threaded	318.230s
2.0 mm x 230 mm Steinmann Pin, trocar/ triangle, fully threaded	381.230s
2.5 mm x 230 mm Steinmann Pin, trocar/triangle, fully threaded	382.230s
3.0 mm x 230 mm Steinmann Pin, trocar/ triangle, fully threaded	383.230s

The Steinmann Pins on this list are the most popular sizes in our range. If you don't see what you need on this list, let us know, we may have it in stock. If not, we will be happy to manufacture to your specification.

**All our pins are:**

- CE marked and manufactured to international standards
- Individually packed in tough Tyvek® backed polyethylene pouches
- Sterile with a shelf life of five years
- Supplied with three traceability labels

The guide wires on this list are the most popular sizes in our range. They are all held in stock ready for rapid shipment to your stores. If you don't see what you need on this list, let us know, we may have it in stock. If not, we will be happy to manufacture wires to your specification, and to hold stock to meet your future requirements.



wires & pins

## guide wires, single trocar with threaded tip

### Description

### Product Code

1.4 mm x 150 mm guide wire, single trocar, threaded tip	375.150s
1.6 mm x 150 mm guide wire, single trocar, threaded tip	348.150s
2.0 mm x 150 mm guide wire, single trocar, threaded tip	359.150s
2.0 mm x 230 mm guide wire, single trocar, threaded tip	359.230s
2.5 mm x 230 mm guide wire, single trocar, threaded tip	353.230s
3.2 mm x 300 mm guide wire, single trocar, threaded tip	376.300s
3.2 mm x 400 mm guide wire, single trocar, threaded tip	376.400s

bone taps

## bone taps



### Description

### Product Code

#### All our sterile packed taps are:

- Individually packed in tough Tyvek® backed polyethylene pouches
- Sterile with a shelf life of five years
- Supplied with three traceability labels per tap

1.50 mm x 50 mm cortical tap, quick fit	947.050s
2.7 mm x 125 mm cortical tap, quick fit	940.125s
3.5 mm x 110 mm cortical tap, quick fit	941.110s
4.5 mm x 125 mm cortical tap, quick fit	942.125s
4.0 mm x 110 mm cancellous tap, quick fit	945.110s



**KellMed**

## packing and sterilisation

All our products are packed in inner and outer pouches of 75gsm tough Tyvek with a clear polyethylene front, heat sealed and fully manufactured to international standards. All products are sterilised with a 25 kGy gamma radiation dose and have a five year shelf life post sterilisation. Each individual pouch carries three self-adhesive traceability labels which are easy and convenient to apply to patient notes where necessary. Drills are supplied in cartons of five and wires in quantities of ten.

## ordering

To place an order please email, phone or fax with the relevant catalogue reference numbers, quantities required and a purchase order number. Thank you.

### Phone & Fax

Tel: +44 (0) 1305 751821  
Fax: +44 (0) 1305 268920  
email: [office@kellmed.co.uk](mailto:office@kellmed.co.uk)

### Sales Office

KellMed Sales Ltd, Retro House, Parkway Farm Business Park  
Poundbury, Dorchester, Dorset DT1 3AR  
[www.kellmed.co.uk](http://www.kellmed.co.uk)





#### Phone & Fax

Tel: +44 (0) 1305 751821  
Fax: +44 (0) 1305 268920  
email: [office@kellmed.co.uk](mailto:office@kellmed.co.uk)

#### Sales Office

KellMed Sales Ltd, Retro House, Parkway Farm Business Park  
Poundbury, Dorchester, Dorset DT1 3AR  
[www.kellmed.co.uk](http://www.kellmed.co.uk)